Mass spectrometry Proteomics and MIAPE

Florian Breitwieser - Research Center of Molecular Medicine, Vienna, Austria

Nov 18, 2010

Outline

BioC-devel meeting Europe 17.-18. 11. 2010 MS-MS Proteomics Standard classes Protein Grouping

Detection of proteins using mass spectrometer

Workflow

- 1. Sample: Proteins
- 2. Digest to: Peptides (Less diversity, more complexity)
- Mass-Spectrometer detects: m/z & Intensity (MS) + Spectrum (MS/MS)
- 4. Search-engine identifies: Peptides and Proteins

Difficulties

- Nature of data
- Standardisation: Formats and minimal information (reproducible research)
- Infrastructure: R and Bioconductor

HUPO PSI and MIAPE

- ► HUPO Proteomics Standards Initiative
 - ► Founded in April 2002
 - http://psidev.info/
 - ► Aims:
 - Create minimum reporting standards
 - ▶ Enable easier transfer of proteomics data
- MIAPE
 - Minimum Information about a proteomics experiment:
 - Mass Spectrometry and Mass Spectrometry Informatics
 - Gel Electrophoresis
 - Liquid Chromatography
 - ► Formats for MIAPE:MS (2.24) / MIAPE:MSI (1.1):
 - mzML for spectrum level information
 - mzldentML for protein identification

Proteomics Formats

- Spectrum formats:
 - ▶ dta, pkl, mgf

 - Seattle Proteome Center: mzXML
 - Proteomics Standard Initiative: mzData,mzML
- ► Identification formats
 - ► Mascot: DAT files
 - Phenyx: PIDRES XML files
 - X Tandem: XML files

 - SPC: pepXML, protXML
 - ▶ PSI: mzldentML

Classes: MIAPE:MS Experiment Information

```
setClass("MIAPE".
         representation=representation(
           creationDate="character".
           contact="character".
           software="character".
           software.version="character".
           software.contact="character".
           software.constomisations="character",
           software.uri="character".
           spectra.source.files="character",
           spectra.source.format="character",
           search.database="character",
           search.database.version="character",
           search.database.nSeg="character".
           search.database.filters="character".
           search.database.nSegSearched="character",
           search.enzvmes="character".
           search.enzymes.missedcleavages="numeric",
           search.enzymes.aditionalParams="character",
           search.modif.fixed="character".
           search.modif.variable="character",
           search.param.fragmentTolerance="character",
           search.param.parentTolerance="character",
           search.threshold.protein="list",
           search.threshold.peptide="list",
           search.additionalParams="list")
```

Classes: MIAPE:MS ProteinGroup and Mzldent

```
setClass("ProteinGroup",
         representation(
                        proteinDescription = "data.frame",
                        ##AC.Description, Validation Status, #different peptide seq, coverage%
                        peptideToProtein = "matrix",
                        peptideDetails = "data.frame",
                        indistinguishableProteins = "character".
                        proteinGroups = "data.frame"
setClass("MzIdent",
         contains = "eSet",
         representation(
           proteinGroup = "ProteinGroup",
           assayData = "list",
                                                # 454+
           ## spectra-matrices: columns m/z, intensity, charge
           phenoData = "AnnotatedDataFrame". # eSet
           featureData = "AnnotatedDataFrame", # eSet
           ## retention time, peptide sequence, peptide modif, scores,
           ## charge state, calculatedMassToCharge, experimentalMassToCharge, ...
           experimentData = "MIAPE",
           annotation = "character",
                                                # eSet
           protocolData="AnnotatedDataFrame"
                                               # eSet
```

Protein Grouping

Problem

- Peptides are detected, not proteins
- Peptides are often shared between proteins
 - especially between splice variants!
- Create a minimal set: Occams Razor

Algorithm

- For each protein: assemble list of peptides with which it is identified
- 2. Indistinguishable proteins: Proteins detected with the same peptides
- 3. Master proteins: Proteins with peptides specific to them
- 4. Group proteins with no specific peptides to master proteins
- 5. Classify

Proteomics in Bioconductor

- ► Good class representations
 - MIBBI: Minimum Information for Biological and Biomedical Investigations
- Multiple experiments?
- Quantitation
- ► Importers:
 - mzML and mzldentML
 - ▶ pepXML, mgf, . . .
- -> Statistics
- -> Applications